

60 85 100 105

	Real 1	Real 2	Real 3	
Speech Understanding	MoreSound Intelligence™ 2.0	Level 1	Level 2	Level 3
	- Environment configuration	5 Options	5 Options	3 Options
	- Virtual Outer Ear	3 Configurations	1 Configuration	1 Configuration
	- Spatial Balancer	100%	60%	60%
	- Neural Noise Suppression, Difficult/Easy	10 dB/4 dB	6 dB/2 dB	6 dB/0 dB
	- Sound Enhancer	3 Configurations	2 Configurations	1 Configuration
	- Wind & Handling Stabilizer	•	•	•
	MoreSound Amplifier™ 2.0	•	•	•
	- SuddenSound Stabilizer	6 Configurations	5 Configurations	4 Configurations
	Feedback Prevention	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield
	Spatial Sound™	4 Estimators	2 Estimators	2 Estimators
	Soft Speech Booster	•	•	•
Sound Quality	Frequency lowering	Speech Rescue™	Speech Rescue™	Speech Rescue™
	Clear Dynamics	•	•	-
	Better-Ear Priority	•	•	-
	Fitting Bandwidth*	10 kHz	8 kHz	8 kHz
	Bass Boost (streaming)	•	•	•
	Processing Channels	64	48	48
	Fitting Bands	24	20	18
Personalisation & Optimising Fitting	Multiple Directionality options	•	•	•
	Adaptation Manager	•	•	•
	Fitting Formulas	VAC+, NAL-NL1/NAL-NL2, DSL v5	VAC+, NAL-NL1/NAL-NL2, DSL v5	VAC+, NAL-NL1/NAL-NL2, DSL v5
Connecting to the world	Oticon Companion app	•	•	•
	Hands-free communication**	•	•	•
	Direct streaming***	•	•	•
	ConnectClip	•	•	•
	EduMic	•	•	•
	Remote Control 3.0	•	•	•
	TV Adapter 3.0	•	•	•
	Phone Adapter 2.0 (with ConnectClip)	•	•	•
	Tinnitus SoundSupport™	•	•	•
	CROS/BiCROS support	•	•	•

* Bandwidth accessible for gain adjustments during fitting

** Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS 15.2 or later

*** From iPhone, iPad, iPod touch, and select Android devices with the Audio Streaming for Hearing Aids (ASHA) protocol