

# PERFORMANCE & FEATURES

	WIDEX MOMENT™ 440	WIDEX MOMENT™ 330	WIDEX MOMENT™ 220	WIDEX MOMENT™ 110
<b>AUDIBILITY/LISTENING IN QUIET</b>				
Widex PureSound™ with ZeroDelay™ Technology	●	●	●	●
TruAcoustics™ with Personal Gain Integrator	●	●	●	●
Variable Speed Compression	● (IE)	● (IE)	●	●
Digital Pinna	●	●	—	—
Enhanced Dynamic Range Compression	●	●	●	●
Fitting rationale options	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL
<b>LISTENING IN NOISE</b>				
True-Input Technology	●	●	●	●
High Definition Locator	15 ch	12 ch	10 ch	6 ch
Speech Enhancer	● (RT/IE)	● (IE)	—	—
Noise Reduction	●	●	●	●
Soft Level Noise Reduction	●	●	●	●
Directional Focus	Via MOMENT App	Via MOMENT App	Via MOMENT App	Via MOMENT App
Reverse focus program	Via MOMENT App	Via MOMENT App	Via MOMENT App	—
<b>USER COMFORT</b>				
Preference Control	●	●	●	●
Fluid Sound Analyser (sound classes)	11 (IE)	7 (IE)	4	3
Multi-directional Active Feedback Cancelling	IE	IE	IE	●
TruSound Softener	●	●	●	—
TruSound AOC	●	●	●	●
Hammer Detector	●	●	●	●
SMARTWIND MANAGER™/Wind Noise Attenuation	●	—	—	—





	WIDEX MOMENT™ 440	WIDEX MOMENT™ 330	WIDEX MOMENT™ 220	WIDEX MOMENT™ 110
<b>USER CONTROLS</b>				
SoundSense Adapt	●	●	●	—
Listening programs	5	4	3	3
InterEar volume/program shift	●	●	●	●
Speech (S)/Music (M) streaming template	S/M	S/M	S/M	—
Phone+	●	●	—	—
SmartSpeak/SmartTone	●	●	●	●
Personal Zen	●	●	●	●
Zen/Zen+	IE	IE	IE	IE
Wireless CROS compatible	●	●	●	●
<b>HEARING HIGH FREQUENCIES</b>				
High-Frequency Boost	●	—	—	—
Audibility Extender	●	●	●	●
<b>FITTING AND FINE TUNING</b>				
Processing channels	15	12	10	6
Fine-tuning channels	15	12	10	6
High-resolution Sensogram (9 fitting bands)	●	●	●	●
Sound Diary/Datalogging	●	●	●	●
<b>APPS, ACCESSORIES AND REMOTE CONTROLS</b>				
WIDEX MOMENT™ App with SoundSense Learn*	●	●	●	●
WIDEX TONELINK™ App	●	●	●	●
COM-DEX App	●	●	●	●
TV PLAY™**	●	●	●	●
COM-DEX / COM-DEX REMOTE-MIC	●	●	●	●
RC-DEX	●	●	●	●
TV-DEX	●	●	●	●
FM+DEX	●	●	●	●

\* mRIC R D and RIC 312 D only, no SoundSense Learn in 110

\*\* Available with WIDEX MOMENT™ compatibility late 2020.



# WIDEX MOMENT™ COMPARISON CHART

WIDEX MOMENT™ hearing instruments are available in a variety of technology levels, each designed to accommodate your specific listening requirements.


























	Natural Sound Experience	Speech Perception	Wind Noise Reduction / Outside	Automatic Sound Adjustment	Hearing in Noise (Parties / Social)	Listening to Music	Using the Phone	Car / Transport	Personal Sound	Varied Environment	High Frequencies
WIDEX MOMENT™ 440	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
WIDEX MOMENT™ 330	★★★★	★★★★	★★	★★★★	★★★★	★★★	★★★★★	★★★★	★★★★	★★★★	★★
WIDEX MOMENT™ 220	★★★	★★★	★★	★★★	★★★	★★★	★★★	★★★★	★★★	★★★	★★
WIDEX MOMENT™ 110	★★	★★	★★	★★	★★	★★★	★★★	★★	★★	★★	★★

WIDEX MOMENT™ hearing aids are packed full of the latest technology to provide the most natural sound ever heard in a hearing aid.

	BENEFIT / USE	440	330	220	110
SMARTWIND MANAGER™	Makes speech and other sounds audible in wind - enjoy conversations in the great outdoors.	●	-	-	-
High Frequency Boost	Clear and crisp sound in high frequencies (in music, for example).	●	-	-	-
Digital Pinna	Mimics the normal pinna effect to help locate where a voice is coming from and hear from behind.	●	●	-	-
Speech Enhancer	Reduces background noise and enhances speech - follow the conversation even in noisy environments.	●	●	-	-
Phone+	Cancels the mic on your opposite ear when using the phone to reduce background noise. Great for landline phones or when streaming is unavailable.	●	●	-	-
Soundsense Learn / Adapt	Widex MOMENT™ App uses artificial intelligence to guide you to optimal, personalised hearing. Remembers your preferences for better sound in every listening situation.	●	●	●	-
HD Locator	Better speech understanding in noise.	15ch	12ch	10ch	-
Sound Classes	Automatically adapts to different listening environments, for optimal clarity, comfort and audibility.	11	7	4	3
Programs	Select the programs you would like your hearing aid to have and manually adjust to appropriate program for different environments.	5	4	2	2

# WIDEX MOMENT™ SOUND CLASS TECHNOLOGY QUICK GUIDE

Exclusive Sound Class Technology allows **WIDEX MOMENT™** hearing aids to adapt automatically to different situations - delivering optimal clarity, comfort and audibility in every environment.

SOUND CLASS	LISTENING ENVIRONMENT	USE / BENEFIT	440	330	220	110
MUSIC	Listening to music	Enjoy your favourite music as it's meant to be heard.				
MUSIC: Classical or contemporary	Optimises sound for different types of music	Settings alter for different levels of intensity and variations in sound level / frequency, associated with different music genres.				
QUIET	Quiet environments with limited background noise (e.g. library)	Hear softer background noises for a more natural hearing experience.				
QUIET with Speech	Optimises speech sounds in quiet environments	Hear conversations in quiet environments comfortably, even where speech is soft or hushed.				
PARTY	Where there is lots of background noise / conversation (e.g. Party)	Distracting background noises are less audible.				
PARTY with Speech	Optimises speech sounds in noisy social situations	Focus on the conversation you want to hear in a noisy social situation.				
TRANSPORT	Situations with low-frequency sounds (e.g. car noise)	Reduces audibility of continuous low-frequency sounds.				
TRANSPORT with speech	Optimises speech sounds in situations with low-frequency sounds	Hear the conversation when travelling.				
URBAN	Where there are variable sounds (e.g. supermarket, city, noisy workplace)	Enjoy the sounds of the city naturally and comfortably without becoming disorientated or overwhelmed.				
URBAN with speech	Optimises speech sounds in urban environments	Comfortably hear the conversation without over amplification of variable background noises.				
SOCIAL	Where multiple conversations occur simultaneously (e.g. busy restaurant)	Small social gatherings where multiple people are talking around you. Focus on the conversation, without distraction.				

## WIDEX MOMENT™ KEY FEATURES

