



# EVOKE™ COMPARISON CHART

| PRODUCT   | SOUND STREAMING              |                 | EASY     |        |       |           | LISTENING ENVIRONMENT |       |          |            |                |         |         | VARIED       |
|-----------|------------------------------|-----------------|----------|--------|-------|-----------|-----------------------|-------|----------|------------|----------------|---------|---------|--------------|
|           | PERSONAL ADJUSTMENT OF SOUND | SOUND STREAMING | EVERYDAY | FAMILY | TV    | TELEPHONE | CAR                   | MUSIC | SHOPPING | RESTAURANT | OFFICE/MEETING | THEATRE | OUTSIDE | SOCIAL EVENT |
| EVOKE 440 | ★★★★★                        | ★★★★★           | ★★★★★    | ★★★★★  | ★★★★★ | ★★★★★     | ★★★★★                 | ★★★★★ | ★★★★★    | ★★★★★      | ★★★★★          | ★★★★★   | ★★★★★   | ★★★★★        |
| EVOKE 330 | ★★★★★                        | ★★★★★           | ★★★★     | ★★★★   | ★★★★  | ★★★★      | ★★★★                  | ★★★★  | ★★★★     | ★★★★       | ★★★★           | ★★★★    | ★★★★    | ★★★★         |
| EVOKE 220 | ★★★★★                        | ★★★★★           | ★★★      | ★★★    | ★★★   | ★★★       | ★★★                   | ★★★   | ★★★      | ★★★        | ★★★            | ★★★     | ★★★     | ★★★          |
| EVOKE 110 | ★★★★★                        | ★★★★★           | ★★       | ★★     | ★★    | ★★        | ★★                    | ★★    | ★★       | ★★         | ★★             | ★★      | ★★      | ★★           |

## EVOKE™ MODEL RANGE


### EVOKE BTEs

- 123 Autumn Beige
- 081 Tan Silk
- 071 Cappuccino Brown
- 069 Summer Gold
- 119 Silver Grey
- 121 Titanium Grey
- 118 Tech Black
- 068 Pearl White
- 076 Sporty Red
- 136 Deep Blue



### EVOKE ITEs


- 025 Light Beige
- 028 Clay Brown
- 027 Dark Brown



## EVOKE™ APPS



**EVOKE™ APP** - for the EVOKE™ Fusion 2 smartphone-compatible model.



### GUIDES YOU TO PERSONALISED HEARING

Refine your hearing in real life with the smart A/B sound profile comparator, which learns from your choices to find your perfect settings, for every environment; Take control in the moment with discreet, simple, smart controls to adjust settings, create personal programs and even find your lost hearing aid!



**TONELINK™ APP** - for all EVOKE™ models



### DISCREET CONTROL OF KEY FUNCTIONS

Turn your smartphone into a remote control to discreetly change programs and adjust settings. There's no need to carry separate assistive listening devices or to touch the hearing aids.

# EVOKE™ SOUND CLASS TECHNOLOGY QUICK GUIDE

Widex EVOKE hearing instruments are available in a variety of technology levels, each designed to accommodate your patients' specific listening requirements. Our exclusive **Sound Class Technology** allows Widex hearing aids to adapt automatically to different situations - delivering optimal clarity, comfort and audibility.

|   | EVOKE440<br>11 Available<br>Sound Classes<br>INTEREAR | EVOKE330<br>7 Available<br>Sound Classes<br>INTEREAR | EVOKE220<br>4 Available<br>Sound Classes | EVOKE110<br>3 Available<br>Sound Classes |
|---|---|--|--|--|
| <b>SOCIAL</b>                                   |   | -  | -  | -  |
| <b>MUSIC</b>                                    |   |  |  |  |
| <b>MUSIC<br/>classical and<br/>contemporary</b> |   | -  | -  | -  |
| <b>QUIET<br/>with speech</b>                    |   |  | -  | -  |
| <b>QUIET</b>                                    |   |  |  |  |
| <b>PARTY<br/>with speech</b>                    |   |  | -  | -  |
| <b>PARTY</b>                                    |   |  |  |  |
| <b>TRANSPORT<br/>with speech</b>                |   | -  | -  | -  |
| <b>TRANSPORT</b>                                |   |  |  | -  |
| <b>URBAN<br/>with speech</b>                    |   | -  | -  | -  |
| <b>URBAN</b>                                    |   |  | -  | -  |

# EVOKE™ KEY FEATURES



The Evoke range is packed full of the latest in hearing aid technology, providing a listening experience so superior that both new and previous Widex hearing aid users will notice the difference. Hearing becomes more effortless and natural.

| FEATURE NAME                                | BENEFIT / USE   | TECHNOLOGY LEVEL |        |     |     |
|---|---|------------------|--------|-----|-----|
|   |   | 440              | 330    | 220 | 110 |
| <b>Processing Channels</b>                  | For accurate processing and a more precise sound experience and speech intelligibility                                  | 15               | 12     | 10  | 6   |
| <b>Fine-tuning Channels</b>                 | More options for personalised sound to cater to your preferences and hearing needs                                      | 15               | 12     | 10  | 6   |
| <b>SoundSense Learn</b>                     | Machine learning guides you to optimal, personalised hearing  | ●                | ●      | ●   | -   |
| <b>SoundSense Adapt</b>                     | Adaptive learning remembers your preferred settings for better sound in every listening situation                       | ●                | ●      | ●   | -   |
| <b>Fluid Sound Controller</b>               | Automatically adjusts to the optimal settings for every environment in real time for a natural sound                    | ●                | ●      | ●   | ●   |
| <b>Fluid Sound Analyzer (Sound Classes)</b> | Accurately analyses your sound environment to prioritise sound sources and select the appropriate sound class           | 11 (IE)          | 7 (IE) | 4   | 3   |
| <b>SmartWind Manager</b>                    | Makes speech and other sounds audible in wind - enjoy conversations in the the great outdoors                           | ●                | -      | -   | -   |
| <b>High Frequency Boost</b>                 | Clear and crisp sound in high frequencies (in music, for example)   | ●                | -      | -   | -   |
| <b>Digital Pinna</b>                        | Mimics the normal pinna effect to help locate where a voice is coming from and hear from behind                         | ●                | ●      | -   | -   |
| <b>Speech Enhancer Real Time</b>            | Reduces background noise and enhances speech in real time so you can follow the conversation even in noisy environments | RT/IE            | IE     | -   | -   |
| <b>HD Locator</b>                           | Better speech intelligibility in noise  | ●                | ●      | ●   | -   |
| <b>TruSound Softener</b>                    | Takes the edge off impact/impulse sounds for improved listening comfort   | ●                | ●      | ●   | -   |
| <b>Evoke Smartphone app</b>                 | Discreet control of key functions (all levels) / Guides you to personalised hearing (220 and above)                     | ●                | ●      | ●   | ●   |
| <b>Tonelink Smartphone App</b>              | Easy control of key functions   | ●                | ●      | ●   | ●   |
| <b>WidexLink wireless connectivity</b>      | Fast, echo-free transmission of sound between hearing aids and DEX devices  | ●                | ●      | ●   | ●   |
| <b>Soft-level noise reduction</b>           | Unwanted soft background noises are reduced for greater listening comfort   | ●                | ●      | ●   | ●   |
| <b>Noise Reduction</b>                      | Reduces disruptive background noise - hear more of the conversation and the sounds important to you                     | ●                | ●      | ●   | ●   |
| <b>ZEN IE</b>                               | Uses soothing tones and/or noise to reduce the painful effects of tinnitus  | ●                | ●      | ●   | ●   |
| <b>Audibility Extender</b>                  | Hear the full conversation - improves perception of important high-frequency speech sounds                              | ●                | ●      | ●   | ●   |
| <b>Variable Speed Compression</b>           | Automatic gearing for smooth sound and stable audibility  | ●                | ●      | ●   | ●   |