

	Intent 1	Intent 2	Intent 3	Intent 4	
	MoreSound Intelligence™ 3.0	Level 1	Level 2	Level 3	Level 4
<b>Speech understanding &amp; listening ease</b>	Sensor technology	✓	✓	-	-
	Environment classifier	5 configurations	5 configurations	3 configurations	Not adjustable
	Virtual Outer Ear	3 configurations	2 configurations	1 configuration	1 configuration
	Spatial Balancer	100%	60%	60%	40%
	Neural Noise Suppression	12 dB / 6 dB	10 dB / 4 dB	8 dB / 2 dB	6 dB / 0 dB
	Difficult/Easy				
	Sound Enhancer	3 configurations	2 configurations	1 configuration	1 configuration
	Wind & Handling Stabilizer	✓	✓	✓	✓
	MoreSound Amplifier™ 3.0	✓	✓	✓	✓
	SuddenSound Stabilizer	6 configurations	5 configurations	4 configurations	2 configurations
	MoreSound Optimizer™	✓	✓	✓	✓
	Feedback shield	✓	✓	✓	✓
	Spatial Sound™	4 estimators	4 estimators	4 estimators	-
	Soft Speech Booster	✓	✓	✓	✓
Frequency lowering, Speech Rescue™	✓	✓	✓	✓	
<b>Sound quality</b>	Clear Dynamics	✓	✓	-	-
	Better-Ear Priority	✓	✓	✓	-
	Fitting Bandwidth <sup>1</sup>	10 kHz	8 kHz	8 kHz	8 kHz
	Power bass (streaming)	✓	✓	✓	✓
	Processing channels	64	48	48	48
<b>Personalisation &amp; optimising fitting</b>	Fitting bands	24	20	18	14
	Multiple directionality options	✓	✓	✓	✓
	Adaptation management	✓	✓	✓	✓
	Fitting formulas	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5
<b>Connecting to the world</b>	Oticon Companion app	✓	✓	✓	✓
	LE Audio	✓	✓	✓	✓
	Hands-free communication <sup>2</sup>	✓	✓	✓	✓
	Direct streaming <sup>3</sup>	✓	✓	✓	✓
	ConnectClip	✓	✓	✓	✓
	EduMic	✓	✓	✓	✓
	Remote Control 3.0	✓	✓	✓	✓
	TV Adapter 3.0	✓	✓	✓	✓
	Phone Adapter 2.0 (with ConnectClip)	✓	✓	✓	✓
	Tinnitus SoundSupport™	✓	✓	✓	✓
	CROS/BiCROS support	✓	✓	✓	✓

<sup>1</sup>Bandwidth accessible for gain adjustments during fitting. <sup>2</sup>Hands-free communication is available on select devices. <sup>3</sup>From iPhone, iPad, Mac, and select Android devices.